

DEPARTMENT  
OF CHEMISTRY



## Centre for Advanced Electron Spin Resonance

"a state-of-the-art facility for use by  
biochemists, chemists and physicists"

October 2024

# Introductory Lecture



- User Applications
- Instrumentation
- ESR Personnel
- How to use CAESR
- Magnet Hazards, Chemical,
- Microwave, and Laser Safety
- Sample Preparation
- Data Storage & Processing
- User Resources



Clarendon




ICL building

# User Applications



- Diamagnetic  $\rightarrow$  Paramagnetic
  - 20 Å to 120 Å
  - Synthesize & Characterize labels:
    - Nitroxides, Trityl, Gd(III);
    - Redox/ E-chem to open shell states
- Intrinsic Paramagnets
  - Paramagnetic Catalyst Intermediates
  - Metallo-enzyme Mechanisms
  - Single Molecule Magnetism
- Transient Paramagnets
  - Photo-Activated Transient States
  - 1 ns resolution
  - Electrical &/or Optical Det

## Pulsed EDMR



# Jose Goicoechea group (1/2)




## CW-EPR characterization

$$h\nu = g_{eff} \mu_B B_0$$

$$g_{eff} = \sqrt{g_{\parallel}^2 \cos^2 \theta + g_{\perp}^2 \sin^2 \theta}$$

by experiment,  $g_{\parallel} = 1.993$   $g_{\perp} = 2.043$



$$\gamma(e^-) = -658 * \gamma(^1H)$$



neat solid, room temperature

J. McGrady DFT predicts

$$g_{\parallel} = 2.000$$
$$g_{\perp} = 2.023$$



*J. Am. Chem. Soc.* **2014**, 136 (4), 1210–1213

CAESR Induction • Oct. 2024




# Jose Goicoechea group (2/2)





## What about $^{73}\text{Ge}$ ? $I = 9/2$

J. McGrady DFT predicts 4x ax. hfc 3 MHz  
8x eq. hfc -20 MHz



$$\text{Pr}(k) = n! f^k (1-f)^{n-k} / k!(n-k)!$$

$n$  total number of sites  
 $k$  number of occupied sites  
 $f$   $^{73}\text{Ge}$  abundance, 7.76%



# Fraser Armstrong (emeritus)




## NiFe Hydrogenase

- Verify that NiC at low pH is consistent with the literature.

- HYSORE was measured at 2.5 K

- Photolysis of the hydride, followed by annealing, leads to characteristic NiC hyperfine interaction.







Dark-adapted sample, mixture of NiC & NiL

**NiL**  
30 min. of laser pulses  
1 mJ pulse energy  
550 nm wavelength  
20 Hz repetition rate

**NiC**  
90 min. annealing at 200K

**NiC simulation**  
 $A(^1H)=[18.4 \ -10.8 \ -18] \text{ /MHz}$   
Lit. ENDOR of *R. eutropha*, Brecht, *et al.*, JACS 2003, 125, 13075.

## Action Spectrum




# ESR Instrumentation



## 5 research instruments

<u>name</u>	<u>GHz</u>	<u>band</u>	<u>method</u>	<u>location</u>
EMX <sub>MICRO</sub>	9.1 - 9.9	X	CW	ICL F19
E580	9.1 - 9.9	X	CW / Pulsed	ICL F11
	33 - 35	Q	CW / Pulsed	
E680	9.1 - 9.9	X	CW / Pulsed	ICL F12
	92 - 94	W	CW / Pulsed	
E380	9.1 - 9.9	X	CW / Pulsed	Clarendon 020
Krymov	130	D	CW / Pulsed	Clarendon 020

- CW-EPR only
- Excellent SNR
- 2.5 – 300 K with ESR-900
- 100 - 450 K with N<sub>2</sub> heater
- 77 K with N<sub>2(l)</sub> finger dewar
- Automatic goniometer
- Room interlocked for Class 4/3B lasers
- User scheduling



# E680

# ICL F12




- 9.1-9.9 GHz, X-band & 94 GHz, W-band
- 2.5 – 300 K with CF-935 & ESR-900
- OPO laser: 213-1700 nm, 6-100 mJ, 7 ns pulse length 20 Hz
- API & ProDel automation
- Arbitrary Waveform Generator (AWG)
- Room interlocked for Class 4/3B lasers

CAESR Induction • Oct. 2024



# E580

# ICL F11



- X-band and Q-band
- Arbitrary Waveform Generator (AWG)
- 140 W at Q-band
- ProDel & API - programmable
- 2.5 – 300 K with CF-935 & ESR-900 cryostats
- OPO laser 355, 410-2200 nm; 4-13 mJ; 7 ns length, 20 Hz
- Room interlocked for Class 4/3B lasers

# CAESR personnel



**Professor Christiane Timmel**  
<http://research.chem.ox.ac.uk/christiane-timmel.aspx>



**Professor Arzhang Ardavan**  
<https://www2.physics.ox.ac.uk/contacts/people/ardavan>

## Dr Will Myers



*Office* 72651  
ICL F20a

*Lab* 72665  
ICL F12

*Fax* 75410

[william.myers@chem.ox.ac.uk](mailto:william.myers@chem.ox.ac.uk)

<http://research.chem.ox.ac.uk/william-myers.aspx>

**Laser contact: Dr. Kevin Henbest**

- 1. Run CAESR & Services in Chemistry
- 4. Consultative role to students and staff
- 5. Collaborate with Department and External groups
- 6. Assist in joint publications
- 7. Pursue, publish, and present peer-reviewed research
- 9. Involved in the submission of grant applications
- 14. Students and staff training
- 15. Responsible for facility safety regulations
- -- Assisting w/ supervision of Part II and D. Phil. student projects
- -- Booking Meetings
- Teaching: Lecturing, small-group teaching, tutoring

**CAESR Induction • Oct. 2024**

# How to Use the CAESR SRF



1. **Independent:** Previous ESR (EPR) or support of experienced group member, ask for any help.
2. **Semi-independent:** Complete training, then independent measurements, scheduling, and data analysis. User asks for support as required.
3. **With CAESR staff:** Complete training, then semi-independent measurements and scheduling. Extended support with project development, data analysis and publications.
4. **Service:** “One-off” projects are rare. Routine or repetitive measurements are generally not possible, (see #2).
5. **Collaboration:** D. Phil student of the Timmel or Ardavan groups &/or Will Myers, works on your project, with measurements, data analysis, & publication(s) as a part of their research.

# Typical User Sequence for ESR



- Discussion with Scientific Applications Manager
  - Research plan, safety, and literature
  - Preliminary sample list
  - Purchase of clear-fused Quartz EPR tubes
- Booking CW-EPR time, COSHH
- Sample Preparations in F13
- CW-EPR in F19, with any required training
- Monthly booking meeting for Pulsed EPR
- Pulsed EPR and / or W-band CW-EPR

# Safety in the ESR Labs.



## Strong Magnetic Fields


Hazards to:

- heart pacemakers
- magnetic back or ID cards
- watches (non-LCD)

Stray fields in Corridors

- outside F12 at a max field

Electromagnet field is mostly within its yoke





# Cryogen and Chemical Safety



- Wet lab work in F13 → departmental guidelines
- All samples ( hazardous or not ), need COSHH form on file in F20a.
- Samples occasionally break -- act quickly to save before complete thawing.
- Safety courses are available for how to handle liquid nitrogen and liquid helium.


# LASER Safety



***Not everyone, but next year ...***

- All CAESR labs are inter-locked for use of Class 4 and 3B lasers.
- Pulsed OPO lasers for 230 - 2200 nm, 4 - 100 mJ, 20 Hz available for Transient ESR. CW lasers at 355 and 532 nm are available.
- New laser users must attend Chemistry Department safety lecture (11:15 am, *next*).
- Users who wish to use their own lasers in CAESR must follow Departmental laser installation procedures.
- Laser Safety forms for each lab must be signed.

# Microwave Safety



$\lambda \sim 3.2 \text{ cm}$

- Place microwave sources in Standby when changing resonators and samples.
- DON'T stare into the end of an open waveguide with microwaves propagating outward.
- X-band (9.4 GHz) exposure may result in cataracts.
  - approx. 60 min of CW at 200 mW / cm<sup>2</sup>

# CAESR Lab ICL F13



- Drying oven, balance, sonicator, vortex mixer
- Fume Cupboard and Schlenk/Vacuum Line
- Flame sealing tubes after de-gassing via FREEZE-PUMP-THAW.
- Flame-sealing J. Young tube bottoms.
- Mbraun LABstar ; Cary60 UV-vis

# Sample tubes and volumes



Tube Dimensions			Micro- wave Band	Sample Vol. ( $\mu\text{L}$ )	<u>Comment</u>
O.D. (mm)	I.D. (mm)	L (mm)			
0.84	0.7	100	W	3	General, for powders and frozen solutions
1.6	1.2	100	Q	10	Dielectric ENDOR resonator
3	2.6	150	X	40	For split ring resonator
3.80	2.79	150	X	80	For ENDOR pulsed resonator, high precision wall thickness
4	3	250	X	250	For all CW resonators and MD5, low-dielectric solutions/solids only

*Tubes of are for sale in ICL F20a, and on the R12 Oracle purchasing system.*



# Sample Concentrations




Transition metals      1 - 5 mM

Organic radicals      20 - 200 mM

Solids       $1 \times 10^{15}$  to  $1 \times 10^{18}$  spins

**Sample** concentration vs. **Spin** concentration

Single Crystals –  $1 \times 10^9$  – auto goniometer



Equal concentrations, Equal noise,  
Different ESR signal width

# Data Storage and Transfer



## Storage

**On the Spectrometers:** Data is kept in  
.../xeprfiles/data/YourPI\_group/You/DateName/...

**Back-Up:** \\chem.ox.ac.uk\SRF\ESR\Spectrometer\_DataBackup\...

→ But no guarantees, you are responsible for your data.

## Transfer

File Transfer Protocol (FTP) is available on E580, E680, and EMXmicro

See: \\chem.ox.ac.uk\SRF\ESR\  
NewUser\_IntrosAndFAQ\FTP\_FileTransfer\_conversion.pdf

# Data Processing Software



Windows, Mac and Linux:

**Matlab** - dept. - Data plotting, integration, spin quantitation, simulations

**Python, C, C++, Java, Fortran, MS Excel**, etc. -free -

**EasySpin** – free, w/ Matlab - [www.easyspin.org](http://www.easyspin.org) – simulation software

**Spinach** – free, w/ Matlab - [www.spindynamics.org](http://www.spindynamics.org) - simulation software

**SpinDynamica** – free, w/ Mathematica - [www.spindynamica.org](http://www.spindynamica.org) - simulation software

+ *many others*

Windows:

**SpinCount** - \$ - Prof. M. Hendrich, Carnegie Mellon Univ. (not in CAESR)

Linux:

**XSophe** - \$\$\$\$ - Bruker BioSpin (not in CAESR)

# User Resources



3<sup>rd</sup> year EPR practical

See: \\chem.ox.ac.uk\SRF\ESR\

Video lecture links

Lecture course slides

Spectrometer Manuals

Stop by -or- questions by e-mail

<http://caesr-web.chem.ox.ac.uk/links.aspx>

- Main Reference Books
  - Pulse EPR methods
  - Selected Systems
  - Instrumentation
- Conferences
- Organizations
- Databases and Software
- Spectrometer Makers
- List of 169 Reference Books

<http://caesr-web.chem.ox.ac.uk/home>

# Graduate Training Sessions



**Morning ~10-12 pm** <https://forms.office.com/e/63ViYPZgXa>

Explanation of webpage and online resources

Tour of ESR Facility

How to Set up cryogenics

How to use a Spectrometer

**Send an interest e-mail to  
william.myers@chem.ox.ac.uk**

Instrumentation description

Demonstration: CW-EPR & Pulsed EPR of LiPc

Pulsed ESR of BDPA: Electron-Nuclear Double Resonance

**Afternoon ~ 2:30 pm to 5 pm**

Demonstration: Transient ESR *with* LASER excitation

Double Electron-Electron Resonance

**Bruker Pulsed EPR Course – group use of the spectrometer**

CAESR Induction • Oct. 2024